

ONOSONIC

Onosonic White Teeth Expert

PRODUCT CATALOG

The Dentist-recommended Premium Brand



Onosonic White Teeth Expert DX-1 ELECTRIC TOOTHBRUSH

The ONOSONIC DX-1 Electric Toothbrush is the result of expert innovation—designed and tested by European expat dentists in Dubai to deliver a professional standard of oral care. Featuring cutting-edge SONIC TWIN ENGINE technology, it produces up to 41,000 strokes per minute for a deep, effective clean that promotes healthier teeth and gums. Thoughtfully crafted with multiple brushing modes to cater to diverse needs, the DX-1 electric toothbrush ensures exceptional care whether you require a gentle touch or a thorough clean. Its built-in smart timer guides you through the dentist- recommended two-minute brushing routine, while the long-lasting battery and ergonomic design make it perfect for both home and travel. Experience the perfect blend of expertise and technology with a toothbrush designed to bring you a confident, radiant smile.



FEATURES

Twin Engine

Generates 41000 strokes per minute 300x more plaque and stain removal

Five Brushing Modes

White: For white teeth and less staining

Clean: Standard daily cleaning

Polish: polishing crowns, veneers, bridges & braces

Gum care: Improves gum health

Sensitive: Gentle cleaning for sensitive areas

Smart Timer

2 minutes timer with 30 seconds pulse indicator

USB Charging

Safe and durable rechargeable lithium battery
Three weeks autonomy

Waterproof

IPX7 waterproof and anti-slip design

TECHNICAL SPECIFICATIONS

| | |
|--|-----------------------|
| Device Model | DX-1 |
| Power Source | DC 3.7V 700mAh |
| Power Consumption | 3.0W |
| CH-416 Charging Station Power Consumption | 1.5W |
| Charging Station Power Source | DC 5V (USB) |
| CH-209 Travel Charging Case | |
| Working Voltage | DC 5V |
| 300ma Power Rating | 1.5 W |
| Battery Capacity | 700mAh |
| Product Material | ABS/TPE/PC/PP |
| Toothbrush Size | 28x26x253mm |
| Package Size | 260x260x90mm |
| IP Rating | IPX7 |

PACKAGE INCLUDES:

- Twin-Engine Brush
- Adult Brush Head
- Junior Brush Head
- Tongue Cleaner
- USB Charger
- USB Travel Case

NOTE:

- Replace the toothbrush head at least every three months.
- Read product instructions carefully when first use.
- The DX-1 is backed by 1-year warranty

Our manufacturer has passed ISO 9001 International Quality System, MDSAP ISO 13485 International Medical Equipment Quality System certification, CCC / CE / ROHS Certification, and other professional institutions' rigorous audits certification.

Onosonic White Teeth Expert DF-1 WATER FLOSSER

The ONOSONIC DF-1 Water Flosser is an advanced oral care device designed to provide a more effective and thorough clean than traditional flossing. Regular use of the DF-1 can significantly improve gum health, reduce plaque buildup, and help prevent cavities. It is recommended to use the device for just two minutes once a day as part of your oral care routine. The DF-1 is particularly beneficial for individuals with braces, implants, crowns, or bridges, as it effectively cleans hard-to-reach areas around orthodontic appliances. With its powerful water jet technology, the DF-1 ensures a deep, refreshing clean, promoting healthier gums and a fresher mouth every day.



FEATURES

Three Modes

High
Daily
Pulse

Nozzles

Standard: for oral rinse food residue
Tongue scraper: for tongue coating cleaning
Periodontal pocket tip: for periodontitis & sensitive gum
Orthodontic: for braces or dentures

USB Charging

Safe and durable rechargeable lithium battery.
30 times, daily mode (2 minutes)

Waterproof

IPX7 waterproof and anti-slip design

Timer

2-minute auto-off

TECHNICAL SPECIFICATIONS

| | |
|-------------------------|---|
| Device Model | DF-1 |
| Material | ABS (Body), PC (Water Tank), PMMA (Nozzles) |
| Input voltage | 5V |
| Power Consumption | 2W |
| Hours of Continuous Use | 30 times, daily mode (2 minutes) |
| Reservoir Capacity | 150 ml |
| Body Weight | 300g |
| Body Size | 63X83X298 mm |
| Rechargeable Battery | Lithium-ion |
| Charge mode | DC (0-50°C/Under Standard atmospheric pressure) |
| Applicable altitude | Below 5000m |
| Charging time | 6 hours |
| IP Rating | IPX7 |

PACKAGE INCLUDES:

- 5 Nozzles Standard X 2 Tongue scraper X1 Periodontal pocket tip X1 Orthodontic X1
- Main Body
- Charging Cable
- Travel Bag

NOTE:

- Read product instructions carefully when first use.
- The DF-1 is backed by 1-year warranty.

Our manufacturer has passed ISO 9001 International Quality System, MDSAP ISO 13485 International Medical Equipment Quality System certification, CCC / CE / ROHS Certification, and other professional institutions' rigorous audits certification.

Onosonic White Teeth Expert **ONO-F1 WATER FLOSSER**

DESCRIPTION

The ONOSONIC WHITE TEETH EXPERT® ONO-F1 Water Flosser delivers a powerful, precision clean that removes plaque and debris from deep between teeth and below the gumline—where toothbrushes and string floss often miss. Ideal for braces, implants, and sensitive gums, it helps reduce bleeding and inflammation while promoting healthier circulation and stronger gum tissue.

With its cordless design, ergonomic grip, and large 300ml detachable water tank, the ONO-F1 ensures a full-mouth clean in one seamless, efficient session. Quiet, travel-friendly, and easy to use, it transforms daily oral care into a refreshing and effective experience.

FEATURES

Adjustable Pressure Settings

Easily control the water pressure by pressing the mode button on the handle. Start with a gentle setting and increase as needed for a comfortable and effective clean.

Nozzles for Every Need

Standard Nozzle: Great for daily use and removing food between teeth.

Tongue Scraper: Helps clean the tongue and keep breath fresh.

Periodontal Tip: Gentle cleaning for sensitive gums or gum pockets.

Orthodontic Nozzle: Designed to clean around braces, crowns, and implants.



USB-C Charging

Comes with a USB-C port for quick and easy charging anywhere.

Waterproof Design

Rated IPX7 waterproof, safe to use at the sink or in the shower, and easy to rinse after use.

Large 300 ml Water Tank

The spacious tank holds enough water for a full-mouth clean in one session—no need to stop and refill.

Smart Timer

A built-in timer helps guide proper cleaning time and turns off automatically when done.

TECHNICAL SPECIFICATIONS

| | |
|-----------------------------|---------------------------------|
| Device Model | ONO-F1 |
| Reservoir Capacity | 300 ml |
| Water Pressure Range | 40-130 psi |
| Input | 5V-1A (USB-C) |
| Charging Time | 3-5 hours |
| Battery Life | 30 + cycles (2 min each) |
| Waterproof Rating | IPX7 |
| Water Flow Method | Pulsating jet |
| Body Materials | ABS (handle), PC (tank) |
| Net Weight | 300 g |

PACKAGE INCLUDES:

- ONO-F1 main handle
- 4 × Nozzles:
 - Standard Jet
 - Orthodontic Jet
 - Periodontal Pocket Tip
 - Tongue-Scraper Tip
- USB-C charging cable
- Travel Poch
- User manual

NOTE:

Our manufacturer has passed ISO 9001 International Quality System, MDSAP ISO 13485 International Medical Equipment Quality System certification, CCC / CE / ROHS Certification, and other professional institutions' rigorous audits certification.

Onosonic White Teeth Expert BRUSH HEAD D-D MEDIUM PACK 2 PCS (ADULT/JUNIOR)

The ONOSONIC D-D Brush Heads are designed for use with the ONOSONIC DX-1 Electric Toothbrush, offering a precise and effective clean. Featuring a unique contour design, these brush heads easily reach all areas of the mouth, including hard-to-reach wisdom teeth. Available in both adult and kids' sizes, the pack includes two brush heads, each proven safe and effective through rigorous quality control tests. The D-D Brush Heads remove up to 100% more stains than a manual toothbrush, promoting whiter teeth and healthier gums. For optimal hygiene, replace the brush heads when the blue bristles become visibly worn or after recovering from a cold or flu. Elevate your oral care with the ONOSONIC Brush Heads.

FEATURES

Precision Design:

Meticulously crafted bristles reach all areas of teeth and gums, strategically targeting hard-to-reach spaces, including gaps between teeth and along the gumline.

Gentle Cleaning:

Despite high-performance capabilities, soft and flexible bristles ensure a gentle brushing experience, minimizing the risk of gum irritation or enamel damage while effectively removing plaque and debris.

High-Quality Materials:

Durable and long-lasting bristles made from high-quality materials maintain their shape and effectiveness over time, ensuring consistent performance for an extended period.

Replacement Reminder:

ONOSONIC toothbrush heads feature a replacement reminder with bristles that have a coloured strip. As the colour fades over time, it indicates when it's time to replace the brush head, maintaining optimal cleaning performance.



Brush Head D-D Medium Pack 2 Pcs-Adult



Brush Head D-D Medium Pack 2 Pcs-Junior

Onosonic White Teeth Expert

WHITENING TOOTHPASTE

Transform your oral care routine with ONOSONIC Teeth Whitening Toothpaste, expertly formulated to deliver professional-level whitening and long-lasting protection. Featuring PAP (Phthalimidoperoxycaproic acid), it effectively removes stains and brightens teeth without causing sensitivity. The inclusion of hydroxyapatite helps repair and strengthen enamel, while xylitol works to protect against cavities, ensuring a healthier mouth with every use. This toothpaste also contains potassium citrate, which is known to reduce dentin hypersensitivity, offering relief for those with sensitive teeth. With its powerful blend of ingredients, ONOSONIC Toothpaste fights stains, strengthens enamel, and provides lasting freshness, leaving your smile whiter, healthier, and more confident.



FEATURES

PAP+ Whitening – Effectively removes tough stains and whitens teeth with minimal sensitivity.

Hydroxyapatite Repair – Restores enamel, strengthens teeth and protects against decay.

Cavity Defence – Xylitol prevents bacterial growth and enhances oral health.

Refreshing Peppermint Oil – Leaves your mouth feeling fresh and clean.

Gentle Formula – Suitable for daily use without harming gums or enamel.

Compact Packaging – Ideal for both home and travel use.

INGREDIENT BREAKDOWN

| CHEMICAL NAME | INCI NAME | CAS | EINECS | CONTENT |
|---|---|-------------|------------|---------|
| Glycerol | Glycerin | 56-81-5 | 200-289-5 | 38.05% |
| Water | Aqua | 7732-18-5 | 231-791-2 | 15.0% |
| Hydrated Silica | Hydrated Silica | 10279-57-9 | 238-878-4 | 11.0% |
| Xylitol | Xylitol | 87-99-0 | 201-788-0 | 10.0% |
| PAP | Phthalimidoperoxycaproic Acid | 128275-31-0 | 410-850-8- | 6.0% |
| Sorbitol | Sorbitol | 50-70-4 | 200-061-5 | 5.0% |
| Hydroxyapatite | Hydroxyapatite | 1306-06-5 | 215-145-7 | 3.0% |
| Potassium Citrate | Potassium Citrate | 866-84-2 | 212-755-5 | 3.0% |
| Peppermint oil | Mentha Piperita (peppermint) oil | 8006-90-4 | | 2.0% |
| Ammonium Acryloyldimethyl-Itaurate/VP Copolymer | Ammonium Acryloyldimethyl-Itaurate/VP Copolymer | | | 1.5% |
| Polyvinylpyrrolidone | PVP | 9003-39-8 | | 1.0% |
| Sodium Hexamethaphosphate | Sodium Hexamethaphosphate | 10124-56-8 | 233-343-1 | 1.0% |
| Piroctone Olamine | Piroctone Olamine | 68890-66-4 | 272-574-2 | 0.9% |
| Titanium Dioxide | Titanium Dioxide | 13463-67-7 | 236-675-5 | 0.9% |
| Sodium Lauroyl Sarcosinate | Sodium Lauroyl Sarcosinate | 137-16-6 | 205-281-5 | 0.9% |
| Mica | Mica | 12001-26-2 | 310-127-6 | 0.25% |
| Saccharin Sodium | Saccharin Sodium. | 6155-57-3 | 204-886-1 | 0.2% |
| Arginine | Arginine | 74-79-3 | 200-811-1 | 0.2% |
| Potassium Hydroxide | Potassium Hydroxide | 1310-58-3 | 215-181-3 | 0.1% |

ONOSONIC Teeth Whitening products have passed rigorous testing and are clinically proven effective, with approval from the relevant regulatory authorities in Dubai.

Onosonic White Teeth Expert **COLOUR CORRECTOR**

DESCRIPTION

The ONOSONIC Colour Corrector is an innovative oral care solution that harnesses advanced colour correcting technology to instantly brighten your smile. This purple-tinted formula works by applying the principles of colour theory—utilizing purple as the complementary colour to yellow on the colour wheel to neutralize unwanted yellow undertones in your teeth. The water-soluble purple dye effectively conceals stains and enhances brightness, providing immediate visual improvement to your smile. Perfect as a post-whitening treatment, the V34 can also be used daily or before special events to maintain optimal tooth brightness. This gentle yet effective formula offers a convenient way to achieve a whiter, more confident smile without the need for harsh bleaching agents.

FEATURES

Advanced Colour Correcting Technology

Utilizes purple pigments to neutralize yellow undertones and instantly brighten teeth appearance through colour theory principles.

Water-Soluble Purple Dye

Safe, effective formula that temporarily cancels out yellow stains while being gentle on tooth enamel.

Versatile Usage Options

Can be used as a post-whitening treatment, daily maintenance routine, or quick brightening boost before events.



Instant Results

Provides immediate visual improvement to tooth brightness without waiting periods or lengthy treatments.

Gentle Formula

Contains xylitol for cavity prevention and peppermint oil for fresh breath, making it suitable for daily use.

Professional-Grade Ingredients

Formulated with high-quality components including glycerin for moisture retention and hydrated silica for gentle polishing.

| NO. | FULL INGREDIENT NAME | INCI NAME | CAS # | CONCENTRATION IN PERCENT | FUNCTION |
|-----|---------------------------|---------------------------|------------|--------------------------|-------------------|
| 1 | Glycerol | GLYCERIN | 56-81-5 | 33.00% | Humectant |
| 2 | Water | AQUA | 7732-18-5 | 27.70% | Solvent |
| 3 | Sorbitol | SORBITOL | 50-70-4 | 25.00% | Humectant |
| 4 | Hydrated Silica | HYDRATED SILICA | 10279-57-9 | 3.00% | Abrasive |
| 5 | Xylitol | XYLITOL | 87-99-0 | 3.00% | caries prevention |
| 6 | Polysorbate 80 | POLYSORBATE 80 | 9005-65-6 | 3.00% | surfactant |
| 7 | Carboxymethyl | CELLULOSE GUM | 9004-32-4 | 2.50% | thickener |
| 8 | Peppermint oil | MENTHA PIPERITA OIL | 8006-90-4 | 1.00% | aromatic |
| 9 | Phenoxyethanol | PHENOXYETHANOL | 122-99-6 | 0.50% | preservative |
| 10 | Sucralose | SUCRALOSE | 56038-13-2 | 0.50% | sweetening agent |
| 11 | Tetrasodium Pyrophosphate | TETRASODIUM PYROPHOSPHATE | 7722-88-5 | 0.50% | chelating agent |
| 12 | C117200 | ACID RED 33 | 3567-66-6 | 0.10% | Colorant |
| 13 | CI 42090 | CI 42090 | 3844-45-9 | 0.10% | Colorant |
| 14 | Ethylhexylglycerin | ETHYLHEXYLGLYCERIN | 70445-33-9 | 0.10% | preservative |
| | | | | 100.000 | |

ONOSONIC Teeth Whitening products have passed rigorous testing and are clinically proven effective, with approval from the relevant regulatory authorities in Dubai.

Onosonic White Teeth Expert

WHITENING POWDER PAP+

Elevate your teeth whitening routine with ONOSONIC Whitening Powder, a revolutionary, non-abrasive solution crafted for deep cleaning and effective stain removal. Formulated with PAP, this advanced whitening powder offers impressive results while being gentle on your enamel, making it safe for daily use. The inclusion of calcium carbonate provides thorough polishing for a smooth, polished finish, while cyclodextrin helps prevent future staining, ensuring your smile stays bright and fresh. Whether you're at home or on the go, ONOSONIC Whitening Powder delivers a quick, easy, and highly effective way to enhance the natural radiance of your smile with lasting results.



FEATURES

Advanced PAP Whitening – Removes stains and whitens teeth with precision.

Calcium Carbonate Polishing – Smooths and brightens enamel without abrasion.

Cyclodextrin Protection – Prevents new stains from forming.

Hydroxyapatite Support – Strengthens and repairs enamel.

User-Friendly Design – Simple to use with your favourite toothbrush.

Travel-Friendly Size – Convenient for use at home or while traveling.

INGREDIENT BREAKDOWN

| CHEMICAL NAME | INCI NAME | CAS | EINECS | CONTENT |
|---------------------------|--------------------------------|-------------|------------|---------|
| Calcium Carbonate | Calcium Carbonate | 471-34-1 | 471-34-1 | 61.85% |
| Cyclodextrin | Cyclodextrin | 7585-39-9 | 231-493-2 | 12.0% |
| PAP | Phthalimidoperoxyacaproic Acid | 128275-31-0 | 410-850-8- | 12.0% |
| Xylitol | Xylitol | 87-99-0 | 201-788-0 | 3.0% |
| Sodium Bicarbonate | Sodium Bicarbonate | 144-55-8 | 205-633-8 | 2.0% |
| Aroma | Aroma | 61789-32-0 | 263-052-5 | 2.0% |
| Sodium Cocoyl Isethionate | Sodium Cocoyl Isethionate | | | 2.0% |
| PVP | PVP | 9003-39-8 | 231-791-2 | 2.0% |
| Water | Aqua | 7732-18-5 | | 0.5% |
| Sucralose | Sucralose | 56038-13-2 | 259-952-2 | 0.5% |
| Hydroxyapatite | Hydroxyapatite | 1306-06-5 | 215-145-7 | 0.5% |
| Mica (C177019) | Mica | 12001-26-2 | 310-127-6 | 0.5% |
| Tetrasodium Pyrophosphate | Tetrasodium Pyrophosphate | 7722-88-5 | 231-767-1 | 0.3% |
| L-Menthol | Menthol | 89-78-1 | 201-939-0 | 0.3% |
| Ethylhexylglycerin | Ethylhexylglycerin | 70445-33-9 | 408-080-2 | 0.2% |
| Phenoxyethanol | Phenoxyethanol | 122-99-6 | 204-589-7 | 0.15% |
| Potassium Citrate | Potassium Citrate | 866-84-2 | 212-755-5 | 0.1% |
| Monosodium Citrate | Monosodium Citrate | 18996-35-5 | 242-734-6 | 0.1% |

ONOSONIC Teeth Whitening products have passed rigorous testing and are clinically proven effective, with approval from the relevant regulatory authorities in Dubai.

Onosonic White Teeth Expert **WHITENING GEL PEN**

Achieve a brighter, more confident smile wherever you go with the ONOSONIC PAP Whitening Gel Pen. Powered by a HP formula, this innovative gel pen provides professional-grade whitening results in a convenient, portable format. Its precision applicator targets stains with ease, giving you a radiant smile in just minutes. Infused with potassium nitrate, it minimizes sensitivity, making it suitable for those with more sensitive teeth, while hydroxyapatite works to repair and protect enamel. The refreshing menthol infusion leaves your mouth feeling clean and invigorated. Whether you're at work, traveling, or out with friends, the ONOSONIC Whitening Gel Pen ensures you always have the confidence to flash a brilliantly white smile.



FEATURES

Precision Whitening – Effectively removes stains from all teeth surfaces.

Hydrating, Sensitivity-Friendly Formula – minimize irritation and keep teeth hydrated.

Xylitol-Enhanced Cavity Defense – Helps protect enamel and reduce plaque buildup during whitening.

Refreshing Menthol Flavor – Leaves your mouth feeling clean and fresh.

Portable Design – Compact pen for convenient, on-the-go whitening.

Simple Application – Easy brush tip for precise, mess-free usage.

INGREDIENT BREAKDOWN

| Full Ingredient Name | INCI NAME | CAS | Concentration in percent |
|-------------------------------|-------------------|------------|--------------------------|
| Propylene Glycol | Propylene Glycol | 57-55-6 | 32.00% |
| Glycerol | Glycerin | 56-81-5 | 30.00% |
| Water | Aqua | 7732-18-5 | 24.00% |
| Hydrogen Peroxide | Hydrogen Peroxide | 7722-84-1 | 6.00% |
| Carbomer | Carbomer | 9007-17-4 | 3.00% |
| Carboxymethylcellulose Sodium | Cellulose Gum | 9004-32-4 | 2.40% |
| Xylitol | Xylitol | 87-99-0 | 1.00% |
| Polyvinylpyrrolidone | Pvp | 9003-39-8 | 0.60% |
| Menthol | Menthol | 89-78-1 | 0.50% |
| Sodium Hydroxide | Sodium Hydroxide | 1310-73-2 | 0.30% |
| Mica | Mica | 12001-26-2 | 0.20% |
| | | | 100.00% |

ONOSONIC Teeth Whitening products have passed rigorous testing and are clinically proven effective, with approval from the relevant regulatory authorities in Dubai.

Onosonic White Teeth Expert

WHITENING DRY STRIPS

Effortlessly whiten your teeth with the ONOSONIC HP Teeth Whitening Strips. These innovative, mess-free strips are designed to provide professional-grade whitening with the convenience of at-home use. Powered by hydrogen peroxide, the formula effectively removes stains, revealing a brighter smile. Enhanced with glycerin and coconut oil, the strips provide hydration and comfort during application, while menthol ensures a refreshing aftertaste. Their dry design ensures easy application, making them perfect for use anytime, anywhere.



FEATURES

Hydrogen Peroxide Whitening – Effectively removes stains for visible whitening results.

PVP-Based Adhesion – Ensures strips stay securely in place during use.

Hydrating Glycerin – Protects and hydrates gums during whitening.

Coconut Oil Comfort – Promotes gentle application and added nourishment.

Menthol Refreshment – Leaves a clean, refreshing taste after use.

Mess-Free Design – Strips are easy to apply and remove, perfect for travel or daily use.

INGREDIENT BREAKDOWN

| CHEMICAL NAME | INCI NAME | CAS | EINECS | CONTENT |
|------------------------|------------------------------|-----------|-----------|---------|
| PVP(K30) | PVP | 9003-39-8 | | 33.0% |
| Glycerin | GLYCERIN | 56-81-5 | 200-289-5 | 20.0% |
| Hydroxypropylcellulose | HYDROXYPROPYLCELLULOSE | 9004-64-2 | | 15.0% |
| PVP(K90) | PVP | 9003-39-8 | | 12.8% |
| Water | AQUA | 7732-18-5 | 231-791-2 | 9.2% |
| Hydrogen Peroxide | HYDROGEN PEROXIDE | 7722-84-1 | 231-765-0 | 6.0% |
| Coconut Oil | COCOS NUCIFERA (COCONUT) OIL | 8001-31-8 | 232-282-8 | 3.0% |
| Menthol | MENTHOL | 89-78-1 | 201-939-0 | 1.0% |

ONOSONIC Teeth Whitening products have passed rigorous testing and are clinically proven effective, with approval from the relevant regulatory authorities in Dubai.

ONOSONIC

THANK YOU!

CONTACT US

🌐 www.onosonic.com

✉ Info@onosonic.com

📍 **ONOSONIC**
Emirates Towers, Dubai, United Arab Emirates



SCAN NOW!

•   @onosonic